
D-Phenothrin

产品编号: D55152

CAS: 26046-85-5

分子式: C₂₃H₂₆O₃

纯度: 99.26%

InChi Key: InChIKey=SBNFWQZLDJGRLK-RTWAWAEBSA-N

Smiles: CC(C)([C@@H]1C=C(C)C)[C@@H]1C(OCc1cc(Oc2ccccc2)ccc1)=O

外观: 固体粉末

作用通路: Microbiology/virology

作用靶点: Parasite

溶解性: DMSO: 60 mg/ml (171.21 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: D-Phenothrin is used in veterinary medicine to control pests on animals and to protect crops. D-Phenothrin is an orally effective type II synthetic pyrethroid widely used to kill insects, mosquitoes and human lice.