
IPA-3

产品编号: D55167

CAS: 42521-82-4

分子式: C₂₀H₁₄O₂S₂

纯度: 98%

InChi Key: InChIKey=RFAXLXKIACIUDT-UHFFFAOYSA-N

Smiles: OC1=CC=C2C=CC=CC2=C1SSC1=C2C=CC=CC2=CC=C1O

外观: 固体粉末

作用通路: Cytoskeletal Signaling

作用靶点: PAK

溶解性: DMSO:65 mg/mL (185.5 mM);Ethanol:7 mg/mL (19.97 mM);Water:<1 mg/mL

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: IPA-3 is a selective non-ATP competitive Pak1 inhibitor with IC₅₀ of 2.5 μM, no inhibition to group II PAKs (PAKs 4-6).