
Dioctanoylglycol

产品编号: D55178

CAS: 627-86-1

分子式: C₁₈H₃₄O₄

纯度: 98.01%

InChi Key: InChIKey=HTNFLUQQANUSLR-UHFFFAOYSA-N

Smiles: CCCCCCCC(OCCOC(CCCCCC)=O)=O

外观: 固体粉末

作用通路: Others

作用靶点: Others

溶解性: DMSO:<30mg/ml

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Dioctanoylglycol is a diacylglycerol kinase inhibitor.