
Tiadinil

产品编号: D55180

CAS: 223580-51-6

分子式: C₁₁H₁₀ClN₃O₂S

纯度: 98%

InChi Key: InChIKey=VJQYLJSMBWXGDV-UHFFFAOYSA-N

Smiles: Cc1ccc(NC(=O)c2snnc2C)cc1Cl

外观: 固体粉末

作用通路: Microbiology/virology

作用靶点: Antibacterial

溶解性: DMSO: 11 mg/mL (41.09 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Tiadinil, a plant activator of systemic acquired resistance, boosts the production of herbivore-induced plant volatiles.