
Olorinab

产品编号: D55189

CAS: 1268881-20-4

分子式: C₁₈H₂₃N₅O₃

纯度: 97.13%

InChi Key: InChIKey=ACSQLTBPYZSGBA-GMXVVIOVSA-N

Smiles: CC(C)(C)[C@@H](CO)NC(C1=NN(C2=NC=C[N+](O-))=C2)C3=C1C[C@@]4([H])[C@]3([H])C4=O

外观: 固体粉末

作用通路: GPCR/G Protein

作用靶点: Cannabinoid Receptor

溶解性: DMSO: 180mg/mL (503.6mM; Need ultrasonic)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Olorinab (APD 371) is a highly potent, selective and fully efficacious cannabinoid receptor type 2 (CB2) agonist. With an EC₅₀ of 6.2 nM for hCB2.