
FMOC-L-1,2,3,4-TETRAHYDRONORHARMAN-3-CAR

产品编号: D55226

CAS: 204322-23-6

分子式: C₂₇H₂₂N₂O₄

纯度: 99.66%

InChi Key: InChiKey=IHHULIKSMJSELI-VWLOTQADSA-N

Smiles: OC(=O)[C@@H]1Cc2c(CN1C(=O)OCC1c3cccc3-c3cccc13)[nH]c1cccc21

外观: 固体粉末

作用通路: Others

作用靶点: Others

溶解性: DMSO: 55 mg/mL (125.43 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: FMOC-L-1,2,3,4-TETRAHYDRONORHARMAN-3-CAR is an organic compound, a derivative of THN. It is used as a starting material for the synthesis of various compounds, including pharmaceuticals, agrochemicals, and other organic compounds.