
MDK-4025

产品编号: D55238

CAS: 66774-02-5

分子式: C₁₉H₁₅ClO

纯度: 96.73%

InChi Key: InChIKey=KTVAHLGKTSPDOG-UHFFFAOYSA-N

Smiles: OC(c1ccccc1)(c1ccccc1)c1ccccc1Cl

外观: 固体粉末

作用通路: Membrane transporter/Ion channel; Metabolism

作用靶点: Calcium Channel

溶解性: DMSO:10 mM

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: MDK-4025 is an inhibitor of the high voltage-activated (HVA) Ca²⁺ current in pyramidal neurons.