
Telotristat Etiprate

产品编号: D55256

CAS: 1137608-69-5

分子式: C₂₇H₂₆ClF₃N₆O₃·C₉H₉NO₃

纯度: 99.21%

InChi Key: InChIKey=XSFPZBUIBYMVEA-CELUQASASA-N

Smiles: OC(=O)CNC(=O)c1ccccc1.CCOC(=O)[C@@H](N)Cc1ccc(cc1)-c1cc(O[C@H](c2ccc(Cl)cc2-n2ccc(C)n2)C(F)(F)F)nc(N)n1

外观: 固体粉末

作用通路: Metabolism

作用靶点: Hydroxylase

溶解性: DMSO:100 mg/mL (132.59 mM);Ethanol:<1 mg/mL;Water:<1 mg/mL

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Telotristat Etiprate (LX 1606 Hippurate) is an orally bioavailable, tryptophan hydroxylase (TPH) inhibitor with potential antiserotonergic activity.