
NS19504

产品编号: D55273

CAS: 327062-46-4

分子式: C₁₀H₉BrN₂S

纯度: 95%

InChi Key: InChIKey=HGWLTZOMQZIUBW-UHFFFAOYSA-N

Smiles: NC1=NC=C(CC2=CC=C(Br)C=C2)S1

外观: 固体粉末

作用通路: Membrane transporter/Ion channel

作用靶点: Potassium Channel

溶解性: DMSO: 249 mg/mL (925.10 mM; Need ultrasonic)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: NS19504 is an activator of Ca²⁺-activated K⁺ channel with EC₅₀ of 11.0 μM. NS19504 shows relaxing effect on bladder smooth muscle spontaneous phasic contractions[1].