
GSK2033

产品编号: D55275

CAS: 1221277-90-2

分子式: C₂₉H₂₈F₃NO₅S₂

纯度: 99.05%

InChi Key: InChIKey=PSOXOVKYGBTPB-UHFFFAOYSA-N

Smiles: Cc(cc1C)cc(C)c1S(N(Cc1ccc(C(F)(F)F)o1)Cc(cc1)ccc1-c1cc(S(C)(=O)=O)ccc1)(=O)=O

外观: 固体粉末

作用通路: Metabolism

作用靶点: Liver X Receptor

溶解性: DMSO: 30 mg/mL (50.70 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: GSK2033 is an antagonist of LXR (pIC₅₀s: 7 and 7.4 for LXR α or LXR β , respectively).