
L-Phenylalanine

产品编号: D10046

CAS: 63-91-2

保存条件: Store at RT.

产品介绍: 级别: BR,

CAS号: 63-91-2

分子式: C₉H₁₁NO₂

分子量: 165.19

MDL: MFCD00064227

溶解性: 50mg/ml 1M Hcl

纯度: ≥99%

外观/性状: 白色结晶

敏感性: 对光敏感

用途: 用于生化研究。配置培养基。营养学研究。

PRODUCT SPECIFIC PUBLICATIONS

[IF=4.1] Rihui Su. et al. Thermosensitive poly(N-isopropylacrylamide) hydrogel/highly internal phase emulsion porous polymer tube tip solid-phase extraction for the determination of methylimidazoles in beverage. J CHROMATOGR A. 2023 Dec;1712:464476 Other ; . 37924617