
(-)-Benzotetramisole

产品编号: D55525

CAS: 950194-37-3

分子式: C₁₅H₁₂N₂S

纯度: 98.22%

InChi Key: InChIKey=YGCWPCVAVSIFLO-GFCCVEGCSA-N

Smiles: N12C=3C(SC1=N[C@H](C2)C4=CC=CC=C4)=CC=CC3

外观: 固体粉末

作用通路: Others

作用靶点: Others

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: (-)-Benzotetramisole makes an efficient and selective catalyst. (-)-Benzotetramisole is a benzannulated derivative of the anthelmintic compound tetraimidazole 3 and is itself an effective catalyst for conversion.