

---

## Telotristat

产品编号: D55539

CAS: 1033805-28-5

分子式: C<sub>25</sub>H<sub>22</sub>ClF<sub>3</sub>N<sub>6</sub>O<sub>3</sub>

纯度: 99.49%

InChi Key: InChiKey=NCLGDOBQAWBXRA-PGRDOPGGSA-N

Smiles: Cc(cc1)nn1-c1c([C@H](C(F)(F)F)Oc2nc(N)nc(-c3ccc(C[C@@H](C(O)=O)N)cc3)c2)ccc(Cl)c1

外观: 固体粉末

作用通路: Metabolism

作用靶点: Hydroxylase

溶解性: DMSO:27.35 mg/mL, with ultrasonic;H<sub>2</sub>O:insoluble;Ethanol:insoluble

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Telotristat is the active metabolite of LX1606(Telotristat etiprate), which is an orally bioavailable tryptophan hydroxylase (TPH) inhibitor and shows potential antiserotonergic activity.