
CTP disodium dihydrate

产品编号: D55870

CAS: 81012-87-5

分子式: C₉H₁₈N₃Na₂O₁₆P₃

纯度: 100%

InChi Key: InChIKey=KNSYTSUGTLJHPS-UHFFFAOYSA-L

Smiles: Nc1ccn(C2OC(COP(=O)(O)OP(=O)([O-])OP(=O)([O-])O)C(O)C2O)c(=O)n1.O.O.[Na+].[Na+]

外观: 固体粉末

作用通路: Membrane transporter/Ion channel;Neuroscience

作用靶点: P2X Receptor

溶解性: DMSO:Insoluble;Water:10 mM

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: CTP disodium dihydrate is an agonist of P2X₄ purinergic receptor