

---

## GPR183

产品编号: D55873

CAS: 1241280-25-0

分子式: C<sub>19</sub>H<sub>19</sub>F<sub>3</sub>N<sub>2</sub>O<sub>2</sub>

纯度: 98%

InChi Key: InChIKey=SIVJKYRAPQKLIM-UHFFFAOYSA-N

Smiles: O=C(NC1=CC(N2CCOCC2)=CC(F)=C1)CCC3=CC=C(F)C(F)=C3

外观: 固体粉末

作用通路: Others

作用靶点: Others

溶解性: DMSO: 50 mg/mL (137.22 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: GPR183 is a chemotactic receptor known for its role in the maturation of B cells, and the endogenous ligand is the oxysterol 7 $\alpha$ ,25-dihydroxycholesterol (7 $\alpha$ ,25-OHC). GPR183 can be used in research on inflammatory bowel disease (IBD).