
GSK-F1

产品编号: D55876

CAS: 1402345-92-9

分子式: C₂₇H₁₈F₅N₅O₄S

纯度: 95.5%

InChi Key: InChIKey=MSRFVAYVUUHQCN-UHFFFAOYSA-N

Smiles: COc(ncc(-c(cc1)cc2c1N=C(N)N(c1c(C(F)(F)F)cccc1)C2=O)c1)c1S(Nc(ccc(F)c1)c1F)(=O)=O

外观: 固体粉末

作用通路: PI3K/Akt/mTOR signaling

作用靶点: PI4K

溶解性: DMSO: 22.5 mg/mL (37.28 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: GSK-F1 is a novel potent PI4KA inhibitor.