
Ximelagatran

产品编号: D55880

CAS: 192939-46-1

分子式: C₂₄H₃₅N₅O₅

纯度: 98.09%

InChi Key: InChIKey=PMOVWDINZOACJB-RBBKRZOGSA-N

Smiles: CCOC(CN[C@H](C1CCCCC1)C(N(CC1)[C@@H]1C(NCc1ccc(C=NNO)cc1)=O)=O)=O

外观: 固体粉末

作用通路: Proteases/Proteasome

作用靶点: Thrombin

溶解性: DMSO: 55 mg/ml (116.14 mM); Water: Insoluble

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Ximelagatran, a direct thrombin inhibitor, can be taken orally and acts solely by inhibiting the actions of thrombin. Ximelagatran is converted to the active agent melagatran in vivo.