

---

## TAK-615

产品编号: D55884

CAS: 1664335-55-0

分子式: C<sub>25</sub>H<sub>22</sub>FNO<sub>4</sub>

纯度: 98%

InChi Key: InChIKey=LEZVXIZCJXKBJY-UHFFFAOYSA-N

Smiles: O=C(O)C1=CC=C(CN(CC2CC2)C(C3=CC=C(OC4=CC=CC=C4F)C=C3)=O)C=C1

外观: 固体粉末

作用通路: GPCR/G Protein

作用靶点: LPA Receptor

溶解性: DMSO : 120 mg/mL (286.10 mM; Need ultrasonic)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: TAK-615 is commonly used to study pulmonary fibrosis and is a negative allosteric modulator of LPA1 receptors.