
Urapidil

产品编号: D55885

CAS: 34661-75-1

分子式: C₂₀H₂₉N₅O₃

纯度: 100%

InChi Key: InChIKey=ICMGLRUYEQNHPPF-UHFFFAOYSA-N

Smiles: COc1cccc1N1CCN(CCCNc2cc(=O)n(C)c(=O)n2C)CC1

外观: 固体粉末

作用通路: GPCR/G Protein; Neuroscience

作用靶点: 5-HT Receptor; Adrenergic Receptor

溶解性: H₂O: < 0.1 mg/mL (insoluble)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Urapidil, a sympatholytic antihypertensive drug, acts as a 5-HT_{1A} receptor agonist and as an α ₁-adrenoceptor antagonist.