
QX-314 chloride

产品编号: D55888

CAS: 5369-03-9

分子式: C₁₆H₂₇ClN₂O

纯度: 99.84%

InChi Key: InChIKey=LLPPOMUAOGMYQI-UHFFFAOYSA-N

Smiles: O=C(NC1=C(C)C=CC=C1C)[N+](CC)(CC)CC.[Cl-]

外观: 固体粉末

作用通路: Membrane transporter/Ion channel

作用靶点: Sodium Channel

溶解性: DMSO:14.29 mg/mL (47.82 mM; Need ultrasonic)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: QX-314 chloride is a membrane-impermeable permanently charged blocker of sodium channel.