
PRX933 hydrochloride

产品编号: D55897

CAS: 639029-42-8

分子式: C₁₆H₂₂ClN₅O₂

纯度: 97.99%

InChi Key: InChIKey=SAQDAMABRIBIHR-BTQNPOSSSA-N

Smiles: C[C@@H]1CNCCN1C2=NC=CN=C2OCCOC3=CC=CN=C3.Cl

外观: 固体粉末

作用通路: GPCR/G Protein;Neuroscience

作用靶点: 5-HT Receptor

溶解性: DMSO: 50 mg/mL (142.11 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: PRX933 hydrochloride is an agonist of 5-HT_{2C} receptor and can be used in studies about hypertension acute treatment.