
Fenebrutinib

产品编号: D55905

CAS: 1434048-34-6

分子式: C₃₇H₄₄N₈O₄

纯度: 98.84%

InChi Key: InChIKey=WNEODWDFDXWOLU-QHCPKHFHSA-N

Smiles: C[C@H]1CN(CCN1C2=CN=C(C=C2)NC3=CC(=CN(C3=O)C)C4=C(C(=NC=C4)N5CCN6C7=C(CC(C7)(C)C)C=C6C5=O)CO)C8COC8

外观: 固体粉末

作用通路: Angiogenesis; Tyrosine Kinase/Adaptors

作用靶点: BTK

溶解性: DMSO:20 mg/mL (30.09 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Fenebrutinib is a selective and noncovalent Btk inhibitor (K_i: 0.91 nM).