
Ciprofibrate

产品编号: D55935

CAS: 52214-84-3

分子式: C₁₃H₁₄Cl₂O₃

纯度: 99.95%

InChi Key: InChIKey=KPSRODZRAIWAKH-UHFFFAOYSA-N

Smiles: CC(C)(OC1=CC=C(C=C1)C1CC1(Cl)Cl)C(O)=O

外观: 固体粉末

作用通路: DNA Damage/DNA Repair; Metabolism

作用靶点: PPAR

溶解性: DMSO:54 mg/mL (186.7 mM);Ethanol:54 mg/mL (186.7 mM);Water:<1 mg/mL

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Ciprofibrate is a fibrate derivative with antilipidemic activity.