
OAC3

产品编号: D55937

CAS: 182564-41-6

分子式: C₁₅H₁₁FN₂O

纯度: 95%

InChi Key: InChIKey=RQJUZWRCVBOMH-UHFFFAOYSA-N

Smiles: O=C(c(cc1)ccc1F)Nc(cc1)cc2c1[nH]cc2

外观: 固体粉末

作用通路: Others

作用靶点: Others

溶解性: DMSO: 60 mg/ml (235.98 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: OAC3 is an activator of octamer-bound transcription factor 4.