
Quinestrol

产品编号: D55938

CAS: 152-43-2

分子式: C₂₅H₃₂O₂

纯度: 98%

InChi Key: InChIKey=PWZUUYSISTUNDW-VAFBSOEGSA-N

Smiles: C[C@]12CC[C@H]3[C@@H](CCc4cc(OC5CCCC5)ccc34)[C@@H]1CC[C@@]2(O)C#C

外观: 固体粉末

作用通路: Endocrinology/Hormones

作用靶点: Estrogen Receptor/ERR; Estrogen/progestogen Receptor

溶解性: DMSO:70 mg/mL (192 mM);Water:<1 mg/mL

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Quinestrol is a synthetic estrogen, applied in hormone replacement therapy, and occasionally to use in cancer and prostate cancer treatment.