
Histone acetyltransferase p300 Inhibitor 4c

产品编号: D55943

CAS: 2103-99-3

分子式: C₉H₇ClN₂S

纯度: 99.86%

InChi Key: InChIKey=DWGWNNCHJPKZNC-UHFFFAOYSA-N

Smiles: NC1=NC(=CS1)C1=CC=C(Cl)C=C1

外观: 固体粉末

作用通路: Metabolism;Neuroscience

作用靶点: Carbonic Anhydrase;AChE

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Histone acetyltransferase p300 Inhibitor 4c is a 2-aminothiazole derivative. Histone acetyltransferase p300 Inhibitor 4c inhibited hCA I, hCA II, AChE and BChE with Ki values of 0.008 ± 0.001, 0.124 ± 0.017, 0.129 ± 0.030 and 0.083 ± 0.041 μM. 0.017, 0.129 ± 0.030 and 0.083 ± 0.041 μM.