
Tocainide

产品编号: D55951

CAS: 41708-72-9

分子式: C₁₁H₁₆N₂O

纯度: 98.04%

InChi Key: InChIKey=BUJAGSGYPOAWEI-UHFFFAOYSA-N

Smiles: CC(N)C(NC1=C(C)C=CC=C1C)=O

外观: 固体粉末

作用通路: Membrane transporter/Ion channel

作用靶点: Sodium Channel

溶解性: DMSO: 11 mg/mL (57.21 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Tocainide is a blocker of the sodium channel and used for the treatment of tinnitus.