
BMS-191095

产品编号: D55955

CAS: 166095-21-2

分子式: C₂₂H₂₁ClN₄O₂

纯度: 98%

InChi Key: InChIKey=SMIKIPXIDLITMP-LEWJYISDSA-N

Smiles: CC(C)([C@@H]([C@H]1N(Cc2ncc[nH]2)c(cc2)ccc2Cl)O)Oc(cc2)c1cc2C#N

外观: 固体粉末

作用通路: Membrane transporter/Ion channel

作用靶点: Potassium Channel

溶解性: DMSO:90 mg/mL (220.11 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: BMS-191095 is mitochondrial ATP-sensitive potassium (mitoKATP) channel activator.