
Drometrizole Trisiloxane

产品编号: D55957

CAS: 155633-54-8

分子式: C₂₄H₃₉N₃O₃Si₃

纯度: 98%

InChi Key: InChIKey=HUVYTMMDZRHBN-UHFFFAOYSA-N

Smiles: CC(Cc(cc(C)cc1-n2nc(cccc3)c3n2)c1O)C[Si](C)(O[Si](C)(C)C)O[Si](C)(C)C

外观: 固体粉末

作用通路: Others

作用靶点: Others

溶解性: DMSO: 5.02 mg/mL (10 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Drometrizole trisiloxane is a photostable UVA and UVB light filter.