
Lirequinil

产品编号: D55963

CAS: 143943-73-1

分子式: C₂₆H₂₅ClN₂O₃

纯度: 97.86%

InChi Key: InChIKey=CBSWRAUYCIIUEI-FQEVSTJZSA-N

Smiles: C(=O)(C1=C2C=3C(CCN2C(=O)C(=C1)C4=CC=CC=C4)=CC=C(Cl)C3)N5C[C@@H](OCC)CC5

外观: 固体粉末

作用通路: Membrane transporter/Ion channel;Neuroscience

作用靶点: GABA Receptor

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Lirequinil is a small molecule GABAA receptor agonist used to study neurological disorders.