
L-Histidine benzyl ester bistosylate

产品编号: D55969

CAS: 24593-59-7

分子式: C₂₇H₃₁N₃O₈S₂

纯度: 99.33%

InChi Key: InChIKey=FQHRMEOHSXSJBX-LTCKWSDVSA-N

Smiles: Cc(cc1)ccc1S(O)(=O)=O.Cc(cc1)ccc1S(O)(=O)=O.N[C@@H](Cc1c[nH]cn1)C(OCc1ccccc1)=O

外观: 固体粉末

作用通路: Others

作用靶点: Others

溶解性: DMSO : 250 mg/mL (423.96 mM; Need ultrasonic)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: L-Histidine benzyl ester bistosylate may serve as an activator for HutP, an RNA-binding protein.