
Bephenium (hydroxynaphthoate)

产品编号: D55987

CAS: 3818-50-6

分子式: C₂₈H₂₉NO₄

纯度: 100%

InChi Key: InChIKey=PMPQCPQAHTXCDK-UHFFFAOYSA-M

Smiles: Oc1cc2ccccc2cc1C([O-])=O.C[N+](C)(CCOc1ccccc1)Cc1ccccc1

外观: 固体粉末

作用通路: Microbiology/virology

作用靶点: Parasite

溶解性: DMSO:198.4 mM

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Bephenium hydroxynaphthoate is a B-type AChR activator, formerly used in the treatment of hookworm infections and ascariasis.