
Temgicoluril

产品编号: D56003

CAS: 10095-06-4

分子式: C₈H₁₄N₄O₂

纯度: 99.83%

InChi Key: InChiKey=XIUUSFJTJXFNGH-UHFFFAOYSA-N

Smiles: CN(C(C(N1C)N2C)N(C)C2=O)C1=O

外观: 固体粉末

作用通路: Membrane transporter/Ion channel; Neuroscience

作用靶点: GABA Receptor

溶解性: DMSO: 125 mg/mL (630.61 mM; Need ultrasonic)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Temgicoluril acts on GABA Receptor.