
Thianthrene

产品编号: D56006

CAS: 92-85-3

分子式: C₁₂H₈S₂

纯度: 99.77%

InChi Key: InChIKey=GVIJJXMXTUZIOD-UHFFFAOYSA-N

Smiles: c(cc1)cc2c1Sc(cccc1)c1S2

外观: 固体粉末

作用通路: Others

作用靶点: Others

溶解性: DMSO: 11 mg/mL (50.85 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Thianthrene is a sulfur-containing heterocyclic chemical compound that is a derivative of the parent heterocycle called dithiin and is easy to oxidize.