
Terbutaline Sulfate

产品编号: D56041

CAS: 23031-32-5

分子式: $C_{12}H_{19}NO_3 \cdot 1/2H_2SO_4$

纯度: 99.47%

InChi Key: InChiKey=KFVSLSTULZVNPG-UHFFFAOYSA-N

Smiles: OS(O)(=O)=O.CC(C)(C)NCC(O)c1cc(O)cc(O)c1.CC(C)(C)NCC(O)c1cc(O)cc(O)c1

外观: 固体粉末

作用通路: GPCR/G Protein; Neuroscience; Microbiology/virology

作用靶点: Adrenergic Receptor; Antibacterial; Antibiotic

溶解性: H₂O:100 mg/mL (364.54 mM; Need ultrasonic); DMSO:5 mg/mL (18.23 mM; Need ultrasonic)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Terbutaline sulfate is a selective beta-2 adrenergic agonist used as a bronchodilator and tocolytic.