
MA-0204

产品编号: D56044

CAS: 2095128-17-7

分子式: C₂₅H₂₇F₃N₂O₄

纯度: 98.83%

InChi Key: InChiKey=GYNMVDMBFKGCCR-QGZVFWFLSA-N

Smiles: C[C@H](CCCOc1c(Cn2c(-c(cc3)ccc3OC(F)(F)F)ncc2C)cccc1)CC(O)=O

外观: 固体粉末

作用通路: DNA Damage/DNA Repair; Metabolism

作用靶点: PPAR

溶解性: DMSO:40 mg/mL (83.95 mM; Need ultrasonic)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: MA-0204 is a highly selective and orally available peroxisome proliferator-activated receptor δ (PPAR δ) modulator (EC₅₀s: 0.4 nM, 7.9 nM and 10 nM for human, mouse and rat PPAR δ , respectively). It is a potential treatment for Duchene Muscular Dystrophy (DMD).