

---

## ISRIB-A15

产品编号: D56050

CAS: 1628478-12-5

分子式: C<sub>22</sub>H<sub>22</sub>Cl<sub>4</sub>N<sub>2</sub>O<sub>4</sub>

纯度: 98.5%

InChi Key: InChIKey=JNLGSKLUABTOBY-HDJSIYSDSA-N

Smiles: O(CC(N[C@H]1CC[C@H](NC(COC2=CC(Cl)=C(Cl)C=C2)=O)CC1)=O)C3=CC(Cl)=C(Cl)C=C3

外观: 固体粉末

作用通路: Others

作用靶点: Others

溶解性: DMSO: Insoluble; H<sub>2</sub>O: Insoluble

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: ISRIB-A15 is an potent inhibitor of the integrated stress response with an EC<sub>50</sub> of 0.8 nM.