
Nedocromil

产品编号: D56053

CAS: 69049-73-6

分子式: C₁₉H₁₇NO₇

纯度: 95%

InChi Key: InChIKey=RQTOOFIXOKYGAN-UHFFFAOYSA-N

Smiles: CCCC(c(OC(C(O)=O)=C1)c(cc23)C1=O)c2N(CC)C(C(O)=O)=CC3=O

外观: 固体粉末

作用通路: GPCR/G Protein; Immunology/Inflammation; Neuroscience

作用靶点: Histamine Receptor; leukotriene receptor; Prostaglandin Receptor

溶解性: DMSO : 16.67 mg/mL (44.89 mM; Need ultrasonic)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Nedocromil inhibits the action or formation of multiple mediators. Which including histamine, leukotriene C₄ (LTC₄), and prostaglandin D₂ (PGD₂).