

---

## 1H-Indole, 3,3'-(2-thienylmethylene)bis[2-methyl-

产品编号: D56063

CAS: 138801-72-6

分子式: C<sub>23</sub>H<sub>20</sub>N<sub>2</sub>S

纯度: 98.8%

InChi Key: InChiKey=FFCNVXLFKIJLEM-UHFFFAOYSA-N

Smiles: S1C=CC=C1C(C=2C=3C=CC=CC3NC2C)C=4C=5C=CC=CC5NC4C

外观: 固体粉末

作用通路: Microbiology/virology

作用靶点: Antifungal

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: 1H-Indole, 3,3'-(2-thienylmethylene)bis[2-methyl- can be used as a broad-spectrum antimicrobial agent with inhibitory activity against *Aspergillus niger* *Bacillus subtilis*, *Pseudomonas albicans*, *Escherichia coli* and *Staphylococcus aureus*. 1H-Indole, 3,3'-(2-thienylmethylene) bis[2-methyl- showed anti-inflammatory activity in in vitro experiments and can be used to treat bacterial infections.