
Transcainide

产品编号: D56076

CAS: 88296-62-2

分子式: C₂₂H₃₅N₃O₂

纯度: 100%

InChi Key: InChIKey=MFKCGXDCHAFQQZ-RTBURBONSA-N

Smiles: C(NC1=C(C)C=CC=C1C)(=O)C2(N(C)C)CCN(CC2)[C@H]3[C@H](O)CCCC3

外观: 固体粉末

作用通路: Membrane transporter/Ion channel

作用靶点: Sodium Channel

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Transcainide(R 54718), a new lidocaine analog, is an orally active antiarrhythmic agent. Transcainide blocks the open state of BTX-activated sodium channels in bovine heart and rat skeletal muscle.