
Dehydroaripiprazole

产品编号: D56079

CAS: 129722-25-4

分子式: C₂₃H₂₅Cl₂N₃O₂

纯度: 97.05%

InChi Key: InChIKey=CDONPRYEWWPREK-UHFFFAOYSA-N

Smiles: O=C1NC2=C(C=CC(OCCCCN3CCN(C4=CC=CC(Cl)=C4Cl)CC3)=C2)C=C1

外观: 固体粉末

作用通路: GPCR/G Protein; Neuroscience

作用靶点: 5-HT Receptor

溶解性: DMSO: 11 mg/mL (24.65 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Dehydroaripiprazole is the active metabolite of aripiprazole. Aripiprazole is an antipsychotic drug, which is metabolized by CYP3A4 and CYP2D6 and mainly forms dehydroaripiprazole. Dehydroaripiprazole has antipsychotic activity equivalent to aripiprazole.