
ML-031

产品编号: D56080

CAS: 852230-33-2

分子式: C₁₉H₂₀N₂O₃

纯度: 98.29%

InChi Key: InChIKey=JRESCQKIRPOOEM-UHFFFAOYSA-N

Smiles: Cc1cc(C(=O)CN2C(=O)CCC2=O)c(C)n1Cc1ccccc1

外观: 固体粉末

作用通路: GPCR/G Protein

作用靶点: S1P Receptor

溶解性: DMSO: 5 mg/mL; Ethanol: 5 mg/mL

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: ML-031 activated S1P2 with an EC₅₀ value of 1 μM in an S1P reporter assay.