
SR11237

产品编号: D56090

CAS: 146670-40-8

分子式: C₂₄H₂₈O₄

纯度: 98.54%

InChi Key: InChiKey=ZZUKALQMHNSTK-UHFFFAOYSA-N

Smiles: CC(C)(CCC1(C)C)c2c1cc(C1(c(cc3)ccc3C(O)=O)OCCO1)cc2

外观: 固体粉末

作用通路: Metabolism

作用靶点: RAR/RXR

溶解性: DMSO : 4 mg/mL (10.51 mM; ultrasonic and warming and heat to 60°C)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: SR11237 is a Pan retinoid X receptor (RXR) agonist. SR11237 causes RXR/RXR homodimers to form and transactivate a reporter gene containing a RXR-response element.