

---

## Exifone

产品编号: D56094

CAS: 52479-85-3

分子式: C<sub>13</sub>H<sub>10</sub>O<sub>7</sub>

纯度: 99.39%

InChi Key: InChIKey=XEDWWPGWIXPVRQ-UHFFFAOYSA-N

Smiles: Oc(ccc(C(c(cc1O)cc(O)c1O)=O)c1O)c1O

外观: 固体粉末

作用通路: Chromatin/Epigenetic; DNA Damage/DNA Repair

作用靶点: HDAC

溶解性: DMSO: 4 mg/mL (14.38 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Exifone is a mixed, nonessential activator of HDAC1 that is capable of binding to both free and substrate-bound enzyme, resulting in an increased relative maximal rate of HDAC1-catalyzed deacetylation.