

---

## PHGDH-inactive

产品编号: D56110

CAS: 1914971-16-6

分子式: C<sub>17</sub>H<sub>21</sub>N<sub>5</sub>S

纯度: 99.66%

InChi Key: InChIKey=ITLPIAAGIQAGRA-UHFFFAOYSA-N

Smiles: Cc1cc(NC(N(CC2)CCN2c2ccncc2)=S)nc(C)c1

外观: 固体粉末

作用通路: Others

作用靶点: Others

溶解性: DMSO: 80mg/ml (244.31mM; Need ultrasonic)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: PHGDH inactive is inactive against PHGDH and serves as a negative control for NCT-502 and NCT-503 [1].