

---

## S1p receptor agonist 1

产品编号: D56130

CAS: 1514888-56-2

分子式: C<sub>23</sub>H<sub>24</sub>FN<sub>3</sub>O<sub>3</sub>

纯度: 99.07%

InChi Key: InChIKey=YBIFMTGYWXNIRZ-UHFFFAOYSA-N

Smiles: c1c(ccc(c1)CC(C)C)c1nc(no1)c1ccc(c(c1)F)CN1CC(C1)C(=O)O

外观: 固体粉末

作用通路: GPCR/G Protein

作用靶点: LPL Receptor; S1P Receptor

溶解性: DMSO:6.8 mg/mL, Need ultrasonic and warming

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: An S1P receptor agonist.