
BMS-394136

产品编号: D56134

CAS: 343246-73-1

分子式: C₂₄H₂₁Cl₂FN₄O

纯度: 98.25%

InChi Key: InChIKey=LADJRPDOKVCQCG-NZQKXSOJSA-N

Smiles: C(=O)(C=1[C@H](N2C(NC1C)=CC=N2)C3=CC(Cl)=C(Cl)C=C3)N4[C@@H](CCC4)C5=CC=C(F)C=C5

外观: 固体粉末

作用通路: Membrane transporter/Ion channel

作用靶点: Potassium Channel

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: BMS-394136 is a KV1.5 antagonist for the treatment of cardiovascular diseases such as arrhythmias and atrial fibrillation.