
Bezafibrate

产品编号: D56136

CAS: 41859-67-0

分子式: C₁₉H₂₀ClNO₄

纯度: 99.82%

InChi Key: InChiKey=IIBYAHWJQTYFKB-UHFFFAOYSA-N

Smiles: CC(C)(Oc1ccc(CCNC(=O)c2ccc(Cl)cc2)cc1)C(=O)O

外观: 固体粉末

作用通路: DNA Damage/DNA Repair; Metabolism

作用靶点: PPAR

溶解性: DMSO:67 mg/mL (185.2 mM);Ethanol:15 mg/mL (41.5 mM);Water:<1 mg/mL

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Bezafibrate is an antilipemic agent that lowers CHOLESTEROL and TRIGLYCERIDES. It decreases LOW-DENSITY LIPOPROTEINS and increases HIGH-DENSITY LIPOPROTEINS.