
GFT505

产品编号: D56140

CAS: 824932-88-9

分子式: C₂₂H₂₄O₄S

纯度: 98.83%

InChi Key: InChiKey=AFLFKFHDSCQHOL-UHFFFAOYSA-N

Smiles: CC(C)(C(O)=O)Oc1c(C)cc(C=CC(c(cc2)ccc2SC)=O)cc1C

外观: 固体粉末

作用通路: DNA Damage/DNA Repair; Metabolism

作用靶点: PPAR

溶解性: DMSO: 60 mg/ml (156.05 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: GFT505 (Elafibranor) is a PPAR α/δ agonist with EC₅₀s of 45 and 175 nM, respectively.