
Iopromide

产品编号: D56145

CAS: 73334-07-3

分子式: C₁₈H₂₄I₃N₃O₈

纯度: 99.4%

InChi Key: InChIKey=DGAIEPBNLOQYER-UHFFFAOYSA-N

Smiles: c1(c(c(c(c1I)NC(=O)COC)I)C(=O)NCC(CO)O)I)C(=O)N(CC(CO)O)C

外观: 固体粉末

作用通路: Others

作用靶点: Others

溶解性: DMSO:91 mg/mL(115 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Iopromide is a low osmolar, non-ionic contrast agent for intravascular use.